

## HB120-S Dry bath

### 제품특징

- A wide range of temperature control up to 120°C
- Timer function support is available
- Various types of optional blocks are available
- Overheating protection
- Lid for protection and heat preservation



### 제품사양

Model	HB120-S
Functions	Heating
Temperature range	Room temp.~120°C
Temperature control accuracy	± 0.5°C
Temperature uniformity	± 0.5°C
Max heating rate	5.5°C/min
Timer	1min-99h59min
Screen	LED
Over temperature protection	140°C
Adapter block material	aluminum
Voltage, Frequency	100-120V/220-240V,50Hz/60Hz
Power	160W
Dimensions[W×D×H ]	175 x 290 x 85mm
weight	3kg

### Blocks for drybath

Cat.No.	Capacity	Diameter	Depth
18900218	0.2mL×54	6.4mm	15mm
18900219	0.5mL×40	8.2mm	26.4mm
18900253 (1.5mm)	1.5mL×40	11mm	34.7mm
18900220 (2mm)	2mL×40	11mm	34.7mm
18900221	5/15mL×28	17mm	48mm
18900222	50mL×8	29.8mm	48mm
18900223	96 or 384 well plate	-	-
18900224	0.2mL+0.5mL+1.5/2mL 18+18+18	-	-