

## dTrite Electronic Titrator

### 제품특징

- High quality motor providing excellent precision and accuracy
- Electronic control decreases the repetitive strain injury
- Remote control panel prevents manual disturbance during operation
- Built-in magnetic stirrer support

### 제품사양

Model	dTrite
Volume range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10 $\mu$ L
Volume accuracy	R=0.2% CV=0.07%
Velocity	16 stages
Operating temperature range	10°C-30°C
Quality control	DIN EN ISO 8655
Control type	External control, including the stirrer and the burette program



## dFlow Electronic Bottle-Top Dispenser

### 제품특징

- Motor operation minimizes the grip force
- Electronic control decreases the repetitive strain injury
- The remote control panel makes the operating process without the manual disturbance
- Two Dispensing modes, dispenser and stepper function options

### 제품사양

Model	dFlow
Volume range	0.1ml-99.9mL, Max piston lift is 10mL, resolution is 100 $\mu$ L
Volume accuracy	Dispenser: R=10mL $\pm$ 30 $\mu$ L; CV=10ml $\pm$ 10 $\mu$ L
Velocity	Stepper R=1mL $\pm$ 6 $\mu$ L; CV=1mL $\pm$ 9 $\mu$ L
Operating temperature range	10°C-30°C
Quality control	DIN EN ISO 8655
Control type	External control

