

EcoStir Economical Magnetic Stirrer

EcoStir is a magnetic stirrer customized for liquid mixing, suitable for chemical analysis, biology reagent mixing, regular laboratory experiments and liquid handling.

제품특징

- Maintenance free brushless DC motor
- Rare earth magnets give powerful stirring
- Stepless speed regulation 300~2000rpm
- Ultra silence and vibration-free design.
- PET top provides excellent chemical resistance
- Shaft can be fitted with a PH sensor or sensors for other purposes



제품사양

Model	EcoStir
Dimension of work plate	φ120mm
Work plate material	PET
Motor type	Brushless DC motor
Motor output power	1.2W
Power	7.2W
Voltage	100-120/200-240V, 50/60Hz
Stirring positions	1
Max. stirring quantity [H ₂ O]	1.5L
Max. magnetic bar [length]	40mm
Speed range	300-2000rpm
Speed setting	scale
Protection class	IP21
Dimension [W × D × H]	150 × 145 × 50mm
Weight	0.3kg
Permissible ambient temperature and humidity	5-40°C, 80%RH