

MS-M-S16 16-Position Magnetic Stirrer

제품특징

- A space saving unit with 16 stirring positions and 4 independent speed control channels
- A wide range of speed from 0 to 1100rpm
- High magnetic adhesion prevents stirrer spin bar out of control
- Stainless steel work plate, covered with silicone cushion, provides excellent performance of skid resistance



제품사양

Model	MS-M-S16
Dimension of work plate[W x D]	365×345mm
Work plate material	Stainless steel with silicone coating
Motor type	DC motor
Motor input power	12W
Motor output power	4W
Power	20W
Voltage	100~120/200~240V,50/60Hz
Stirring positions	16 (4 channel)
Max. stirring quantity [H ₂ O]	0.4LX16
Max. magnetic bar [length]	40mm
Speed range	0-1100rpm
Speed display	scale
Protection class	IP21
Dimension [W x D x H]	380×622×65mm
Weight	6kg
Permissible ambient temperature and humidity	5-40°C, 80%RH