

MS-PA LED Digital Magnetic Stirrer

제품특징

- Digital speed control within the range of 100-1500rpm
- Max. stirring quantity of H₂O at 3L
- LED display shows speed
- Nylon+GF housing provides a good performance of chemical resistance



제품사양

Model	MS-PA
Dimension of work plate	φ135mm (5 inch)
Work plate material	ABS
Motor type	DC motor
Motor input power	5W
Motor output power	3W
Power	15W
Voltage	100-120/200-240V, 50/60Hz
Stirring positions	1
Max. stirring quantity [H ₂ O]	3L
Max. magnetic bar [length]	50mm
Speed range	100-1500rpm
Speed display	LED
Protection class	IP42
Dimension [W x D x H]	150×260×80mm
Weight	1.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH