

MX-RL-Pro LCD Digital Rotator

DLAB rotators give gentle but effective mixing of biological samples held in 1.5mL to 50mL tubes, and are used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics, etc.

제품특징


- Electronic adjustment of speed and time
- LCD display speed and time
- Speed range from 10 to 70rpm
- Provides gentle but effective rolling
- A wide variety of accessories are available



제품사양


Model	MX-RL-Pro
Motor type	DC motor
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Timer	1-1199min
Operation mode	Continuous / timer
Voltage	100-240V, 50/60Hz
Power	40W
Dimensions [D × W × H]	220 × 510 × 260mm
Permissible ambient temperature and humidity	5~40°C, 80%RH
Protection class	IP21

Accessories




Disk Accessories Type-1

Used for MX-RD-Pro
for 1.5mL x 60 centrifuge tubes holder
Cat. No. 18900160



Disk Accessories Type-2

Used for MX-RD-Pro
for 15mL x 16 centrifuge tubes holder
Cat. No.18900161



Disk Accessories Type-3

Used for MX-RD-Pro
for 50mL x 8 centrifuge tubes holder
Cat. No.18900162



Disk support rods

4 pcs in a group
use with disk accessories
Cat. No.18900140



Rotisserie Accessory Type-1

Used for MX-RL-Pro,for 1.5mL x 48
centrifuge tubes held horizontally
Cat. No.18900142



Rotisserie Accessory Type-2

Used for MX-RL-Pro,for 15mL x 24
centrifuge tubes held horizontally
Cat. No.18900143



Rotisserie Accessory Type-3

Used for MX-RL-Pro,for 50mL x 24
centrifuge tubes held horizontally
Cat. No.18900144



Rotisserie Accessory Type-4

Used for MX-RL-Pro,for 1.5mL x 32
centrifuge tubes held vertically
Cat. No.18900145



Rotisserie Accessory Type-5

Used for MX-RL-Pro,for 15mL x 16
centrifuge tubes held vertically
Cat. No.18900146



Rotisserie Accessory Type-6

Used for MX-RL-Pro,for 50mL x 16
centrifuge tubes held vertically
Cat. No.18900147



Plastic Clips

Used with MX-RL-E/MX-RD-E
Disk/Rotisserie accessory, for
50mL, 15mL and 1.5mL
centrifuge tubes