

## MX-T6-Pro LCD Digital Tube Roller

DLAB tube rollers provide a gentle and highly efficient rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, such as the prevention of blood coagulation, immune precipitation, etc.

### 제품특징

- LCD displays the speed and time
- Speed range of 10-70rpm
- Six roller design for gentle rolling and easy to clean
- Continuous or timer operation



### 제품사양

Model	MX-T6-Pro
Motor type	DC motor
Rocking motion	Rocking and rolling
Amplitude	24mm
Max. load capacity	4kg
Number of rollers	6
Roller size [length]	280mm
Speed range	10~70rpm
Speed display	LCD
Timer	Yes
Timer display	LCD
Time setting range	1-1199min
Operation mode	Continuous / timed
Voltage	100-240V, 50/60Hz
Power	30W
Dimensions [D × W × H]	450 × 260 × 120mm
Weight	5.1kg
Permissible ambient temperature and humidity	5~40°C, 80%RH
Protection class	IP21