

SK-D1807-S LED Digital 3D shaker

DLAB rocking and 3D shakers provide perfect solution in electrophoresis gel dyeing/dicoloration, sample washing, molecular hybridization, immunoprecipitation, Western blot and cell culture.

제품특징

- Three-dimensional rocking motion with the angle of 7°
- Max. load capacity is 3kg with platform
- Speed range of 10-80rpm
- Double LEDs for independent display of time and speed
- Various of platform accessories



제품사양

Model	SK-D1807-S
Tilt angle	7°
Max. load capacity(with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Speed range	10-80rpm
Speed display	LED
Timer display	LED
Time settings range	1min-19h59min
Voltage	100-240V, 50/60Hz
Power	20W
Weight	2.7kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Dimensions[W × D × H]	270×330×130mm
Protection class	IP21