

SK-O180-E Economical Digital Orbital Shaker

DLAB orbital and linear shakers are ideal for almost any vessel from microcentrifuge tubes through petri dishes and microtitre plates to conical flasks.

제품특징

- Max. load capacity is 3kg with platform
- Speed range of 0-200rpm
- Wide range of platforms for use with a variety of vessels



제품사양

| | |
|--|------------------|
| Model | SK-O180-E |
| Orbital diameter | 20mm |
| Max. load capacity (with platform) | 3kg |
| Motor type | DC motor |
| Motor input power | 16W |
| Motor output power | 10W |
| Speed range | 0~200rpm |
| Speed display | scale |
| Timer display | - |
| Time settings range | - |
| Voltage | 100~240V,50/60Hz |
| Power | 20W |
| Weight | 3.1kg |
| Permissible ambient temperature and humidity | 5~40°C,80% |
| Dimension[W×D×H] | 270×330×110mm |
| Protection class | IP21 |
| Data connector | - |