

SK-R330-Pro LCD Digital Rocking Shaker

DLAB rocking and 3D shakers provide perfect solution in electrophoresis gel dyeing/dicoloration, sample washing, molecular hybridization, immunoprecipitation, Western blot and cell culture.

제품특징

- Rocking motion angle 9°
- Max. load capacity is 10kg with platform
- Speed range of 10-70rpm
- Double LCDs for independent display of time and speed
- Wide range of platforms for use with a variety of vessels



제품사양

Model	SK-R330-Pro
Tilt angle	9°
Max. load capacity (with platform)	10kg
Motor type	DC motor
Motor input power	40W
Motor output power	24W
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Time settings range	1min-19h59min
Voltage	100-240V, 50/60Hz
Power	50W
Weight	9kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Dimensions [W × D × H]	360 × 410 × 200mm
Protection class	IP21