

T-8 / T-25 ULTRA-TURRAX® 파쇄.균질.분산.유화기



▶ T8 basic



▶ T25 basic

제품특징

Competitively priced dispersing instrument for volumes of 0.5 to 2,000ml (H₂O).
 A wide speed range allows you to work at high circumferential speed.

- Electronic speed control
- Electronic overload protection
- Quick release button for dispersing element

제품사양

Model	T8 basic	T25 basic
Ident No.	8002000	2953000
Motor capacity	100W	500W
Speed range	5,000~25,000 rpm	6,500~24,000 rpm
Speed display	Scale	Scale
Speed adjustment	Stepless	Stepless
Volume range (H ₂ O)	0.5 ~ 50 ml	1 ~ 2,000 ml
Max. viscosity	5,000 mPas	5,000 mPas
Noise without dispersing element	64 dB	73 dB
Overload protection	Yes	Yes
Smooth run-on/start	Yes	Yes
Dimension [W x D x H]	42 x 44 x 250mm	65 x 80 x 240mm
Weight	0.4 kg	1.6 kg